

Dental Department

TrueFace

Manufacturer & Exporter of Surgical &
Dental Instruments

Dental Self Aspirating Syringe

Durable and Best Quality for Perfect Performance.
High Degree of Aesthetic & Corrosion Resistance.

Dental Self Aspirating Anesthetic Syringe. (Thumb Ring)
Available in 1.8ML Injection and 2.2ML Injection



Description;

Made of Quality German Stainless Steel.
Dental Self Anesthetic Syringe.
Best Quality Silver Polish Finish.
for Durability, Perfect & Smooth Working.
Ultrasonic Cleaned.
Autoclave able & Sterilize able.
Pack: Comes in PU Leather Zipper Pouch.



We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument

Dental Department

TrueFace

Manufacturer & Exporter of Surgical &
Dental Instruments

Dental Automatic Syringe

Durable and Best Quality for Perfect Performance.
High Degree of Aesthetic & Corrosion Resistance.



DENTAL INTRALEGAMENTAL SYRINGE

1.8ML Injection. American+English Threads
Precise 0.2cc of Anesthetic with each Click

Description;

Made of Quality German Stainless Steel.
Dental Anesthetic Intralegamental Syringe.
Best Quality Black Colour Finish Handle.
for Durability, Perfect & Smooth Working.
Ultrasonic Cleaned.
Autoclave able & Sterilize able.
Pack: Comes in PU Leather Zipper Pouch.

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Dental Department

TrueFace

Manufacturer & Exporter of Surgical &
Dental Instruments

Dental Automatic Syringe

Durable and Best Quality for Perfect Performance.
High Degree of Aesthetic & Corrosion Resistance.



DENTAL INTRALEGAMENTAL SYRINGE

1.8ML Injection. American+English Threads
Precise 0.2cc of Anesthetic with each Click

Description;

Made of Quality German Stainless Steel.
Dental Anesthetic Intralegamental Syringe.
Best Quality Silver Polish Finish Handle.
for Durability, Perfect & Smooth Working.
Ultrasonic Cleaned.
Autoclave able & Sterilize able.
Pack: Comes in PU Leather Zipper Pouch.

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Dental Department

Dental Drill Guide & Implant Positioning

Made of German Stainless Steel with Best Plastic Box.

TrueFace

Manufacturer & Exporter of Surgical &
Dental Instruments



Dental Drill Guide & Implant Positioning Kit.

Pilot Drill Diameter 2.0mm.

Drill Guide & Guide Pin
Diameters Mentioned Bellow,
6.0 mm
7.5 mm
9.0 mm
10.5 mm
12.0 mm

Best Quality Titanium Anodized Finish for ease in Identification. Passivated, Rusting Prevention Procedure. Ultrasonic Cleaned. Combination of Machine & Hand Crafted, Premium Quality. High Degree of Aesthetic and Corrosion Resistance. Autoclave able & Sterilize able.

We Do Customization

Always follow the Manufacturer's Instructions for Cleaning, Disinfection & Sterilization of Instrument



Dental Department

Dental Drill Guide & Implant Positioning

Made of German Stainless Steel & Titanium Material.

TrueFace

Manufacturer & Exporter of Surgical & Dental Instruments



Dental Drill Guide & Implant Positioning Kit.

Pilot Drill Diameter 2.0mm.

Drill Guide & Guide Pin Diameters Mentioned Bellow,
6.0 mm
7.5 mm
9.0 mm
10.5 mm
12.0 mm

Best Quality Titanium Anodized Finish for ease in Identification. Passivated, Rusting Prevention Procedure. Ultrasonic Cleaned. Combination of Machine & Hand Crafted, Premium Quality. High Degree of Aesthetic and Corrosion Resistance. Autoclave able & Sterilize able.



We Do Customization

Always follow the Manufacturer's Instructions for Cleaning, Disinfection & Sterilization of Instrument

Dental Department

Trephine Drills Kit

Made of German Stainless Steel

Durable and Best Quality for Perfect Performance.
High Degree of Aesthetic & Corrosion Resistance.

TrueFace

Manufacturer & Exporter of Surgical &
Dental Instruments

Dental Implant Trephine Drills Kit.
for Professional Implant Surgery.
Polished Finish, Depth Markings



Trephine Sizes;
(1mm Wall)
1x3.0mm
1x4.0mm
1x5.0mm
1x6.0mm
1x7.0mm
1x8.0mm
1x9.0mm
1x10.0mm

Best Quality Silver Finish.
Passivated.
Rust Prevention Procedure.
Ultrasonic Cleaned.
Combination of
Machine & Hand Crafted.
Autoclave able/Sterilize able.

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Trephine Drills & Saw Disks with Stainless Steel Bur Holder.

Durable and Best Quality for Perfect Performance.
High Degree of Aesthetic & Corrosion Resistance.



Description;
Trephine Sizes:
1 x 3.0mm
1 x 4.0mm
1 x 5.0mm
1 x 6.0mm
1 x 7.0mm
1 x 8.0mm

Saw Disk Sizes;
1 x 7mm
1 x 10mm

Ultrasonic Cleaned.
Autoclave able.
Sterilize able.

Made of German Stainless Steel.

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Bone Expander Saw Disk Kit

with Stainless Steel Bur Holder.

Made of German Stainless Steel.

Durable & Best Quality for Perfect Performance.

High Degree of Aesthetic and Corrosion Resistance.



Description;

Saw Disk Sizes;

1 x 4.0mm

1 x 7.0mm

1 x 10mm

1 x 13mm

Ultrasonic Cleaned.

Autoclave able.

Sterilize able.

Depth Markings.

Precise Cutting: The fine serrated edges allow for clean & controlled cutting of the cortical bone.

Varying Diameters: The Different size diameters allow the surgeon to adapt the instrument to different surgical sites and bone widths along the arch and ensuring a precise cut relative to the anatomy and planned implant size.



We Do Customization

Always follow the Manufacturer's Instructions for Cleaning, Disinfection & Sterilization of Instrument

Dental Department

TrueFace

Manufacturer & Exporter of Surgical &
Dental Instruments

Long Trefphine Drills Kit with Stainless Steel Bur Holder.

Durable & Best Quality for Perfect Performance.
High Degree of Aesthetic & Corrosion Resistance.



Description;
Trimming Burs
Sizes:

- 1 x 3.0mm
- 1 x 4.0mm
- 1 x 5.0mm
- 1 x 6.0mm
- 1 x 7.0mm
- 1 x 8.0mm

Made of German Stainless Steel.
Ultrasonic Cleaned.
Autoclave able and Sterilize able.

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Lateral Approach Drills Kit.

Lateral Approach Sinus Lift Reamer Bur Bone Shape Kit.
with Stainless Steel Drills Holder, Autoclavable & Sterilizable.



Description;

Drill Ø 3.3mm
Drill Ø 3.7mm
Drill Ø 3.3mm
(Head Size)

Drill Ø 8.0mm
Drill Ø 6.0mm
Cutter Ø 8.0mm

Made of German Stainless Steel.
Best Quality Diamond & Titanium Coated for Perfect Performance.
Minimal Invasiveness. Ultrasonic Cleaned.

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Dental Department

TrueFace

Manufacturer & Exporter of Surgical &
Dental Instruments

Universal Drills & Lance Pilot Kit with Stainless Steel Bur Holder.

Durable and Best Quality for Perfect Performance.
High Degree of Aesthetic & Corrosion Resistance.



Description;
Lance Pilot Size:
1 x 2.0mm

Drills Sizes;
1 x 2mm
1 x 2.5mm
1 x 2.8mm
1 x 3mm

Ultrasonic Cleaned.
Autoclave able.
Sterilize able.

Made of German Stainless Steel. **We Do Customization**
Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Universal Drills & Lance Pilot Kit

with Stainless Steel Bur Holder.

Made of German Stainless Steel.

Durable & Best Quality for
Perfect Performance.

High Degree of Aesthetic
and Corrosion Resistance.



Description;

Lance Pilot Size;
1 x 2.0mm

Drills Sizes;

1 x 2.0mm
1 x 2.5mm
1 x 2.8mm
1 x 3.0mm
1 x 3.15mm
1 x 3.3mm
1 x 3.5mm

Ultrasonic Cleaned.
Autoclave able.
Sterilize able.

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Universal Prosthetic Kit

Dental Universal Implant Prosthetic Kit



Description;

Torque Wrench: 10-70 N.cm
Drivers 0.9 (Long + Short)
Drivers 1.0 (Long + Short)
Drivers 1.2 (Long + Short)
Drivers 1.25 (Long + Short)
Drivers 1.27 (Long + Short)
Drivers 1.4 (Long + Short)
Drivers NOB (Long + Short)
Drivers ITI (Long + Short)

1 Torque Wrench.
16 Pcs Screwdrivers.
Short drivers: 8.5mm
Long drivers: 13.5mm

Made of German Stainless Steel.

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument



Dental Department

Dental Bone Expander Kit

Durable & Best Quality for Perfect Performance.
High Degree of Aesthetic & Corrosion Resistance.

TrueFace

Manufacturer & Exporter of Surgical &
Dental Instruments

Best Quality Dental Bone Expander
Kit Lift Saw Disks, Surgical Implant
Instruments.
Ultrasonic Cleaned.
Autoclave able & Sterilize able.



Description;

Made of German Stainless Steel.

Kit Includes;

- 1 x Ratchet Wrench, 4.0mm Square.
- 5 x Bone Expansion Screws
(2.6mm, 3.0mm, 3.4mm, 3.8mm, 4.3mm)
- 1 x Ratchet Extender.
- 1 x 4.0mm Motor Mount Adapter.
- 1 x Hand Driver.
- 1 x 1.8mm Drill
- 2 x Pins
- 1 x Sterilization Cassette - Box

We Do Customization

Always follow the Manufacturer's Instructions for
Cleaning, Disinfection & Sterilization of Instrument

